

## ASPIRA®3P-sPB


### TECHNICAL INFORMATION

MONOFOCAL  
3-PIECE  
ADD-ON



#### ASPIRA®3P-sPB

MS 714 PB



Type

foldable posterior chamber Add-On IOL  
for sulcus fixation in pseudophakic eyes  
three-piece with modified C-loop haptics

Material

optic: silicone with UV absorber, haptics: high molecular PMMA

Specifics

≈ aberration free <sup>(2)</sup>  
convex-concave with round anterior edge for iris protection

Diopter Range

Sphere in 0.5 D steps

-6.0 to 6.0 D

Est. A-const<sup>(1)</sup>

not applicable (by application principle)

<sup>(1)</sup> For more details please visit [www.humanoptics.com](http://www.humanoptics.com)

<sup>(2)</sup> The word "aberration" as used in this document refers to "spherical aberration"