

ASPIRA®3P-sSA


TECHNICAL INFORMATION

MONOFOCAL
3-PIECE
ASPHERIC
SULCUS



ASPIRA®3P-sSA

MS 614 ASP



Type	foldable posterior chamber IOL for sulcus fixation three-piece with modified C-loop haptics
Material	optic: silicone with UV absorber, haptics: high molecular PMMA
Specifics	aspheric posterior surface (aberration correcting) ⁽²⁾
Diopter Range	<i>for more details on diopter categories, see page 2</i>
Sphere in 0.5 D steps (optic shape)	7.0 to 31.0 D (biconvex)
Sphere in 1.0 D steps (optic shape)	32.0 to 40.0 D (biconvex)
Est. A-const ⁽¹⁾	118.0 (optical biometry)

⁽¹⁾ For more details please visit www.humanoptics.com

⁽²⁾ The word "aberration" as used in this document refers to "spherical aberration"

DIOPTER CATEGORY - ASPIRA®3P-sSA

	Standard	7.0 to 31.0 D	in 0.5 D steps
	Extended	32.0 to 40.0 D	in 1.0 D steps