



DIFFRACTIVA® DIFF-sPB

TECHNICAL INFORMATION	
<p>MULTIFOCAL 3-PIECE ADD-ON</p> 	<p>DIFF-sPB MS 714 PB DIFF</p> 
Type	foldable multifocal posterior chamber Add-On IOL for sulcus fixation in pseudophakic eyes three-piece with modified C-loop haptics
Material	optic: silicone with UV absorber, haptics: high molecular PMMA
Specifics	≈ aberration free ⁽²⁾ diffractive, near addition +3.5 D (at IOL plane) convex-concave with round anterior edge for iris protection
Diopter Range	<i>more details on diopter categories see page 2</i>
Sphere in 0.5 D steps (optic shape)	-6.0 to 6.0 D (convex-concave)
Est. A-const ⁽¹⁾	not applicable (by application principle)

⁽¹⁾ For more details please visit www.humanoptics.com

⁽²⁾ The word „aberration“ as used in this document refers to „spherical aberration“

DIOPTRER CATEGORY - DIFFRACTIVA® DIFF-sPB

Standard	0.0 D	
Special	-6.0 to -0.5 D 0.5 to 6.0 D	in 0.5 D steps in 0.5 D steps