



## MARFAN

TECHNICAL INFORMATION	
MONOFOCAL 3-PIECE SULCUS  	<b>MS 714</b>  
Type	foldable posterior chamber IOL for sulcus fixation to patients with Marfan syndrome three-piece with modified C-loop haptics
Material	optic: silicone with UV-absorber, haptics: high molecular PMMA
Specifics	angulation 10°
Diopter Range	
Sphere in 0.5 D steps (optic shape)	7.0 to 31.0 D (biconvex)
Sphere in 1.0 D steps (optic shape)	32.0 to 36.0 D (biconvex)
Est. A-const <sup>(1)</sup>	118.3 (optical biometry)

<sup>(1)</sup> For more details please visit [www.humanoptics.com](http://www.humanoptics.com)