



SAFELOADER®

2-COMPONENT AUTOLOADING SYSTEM



Leading Technology



Premium Quality



Customized Service



Global Experience

INNOVATIVE: THE SAFELOADER® AUTOLOADING CONTAINER

For a relaxed OR atmosphere

The preloading system replaces the manual IOL-loading process into the cartridge

ESTABLISHED: SEPARATED COMPONENTS

For the highest safety level

Separate sterilization processes for the IOL and the injector

Isolated IOL storage:

First contact of the IOL with the cartridge coating immediately before insertion

FLEXIBLE: COMBINE YOUR OWN

Various combinations

The different incision sizes of the Medical Accuject™ push-type injector can be combined with the various Premium IOLs of the monobloc platform


Available sizes of the Accuject™ push-type injector	1.8	2.0	2.2
---	-----	-----	-----




INTUITIVE. SIMPLE. FAST.




Like turning a key



① Insert



② Turn to the right



③ Remove and close

THE PREMIUM MONOBLOC IOL PLATFORM

- **HD optic achieved by sub-nano resolution technology**

Premium-quality optical surface for perfect and brilliant images

- **MICS**

Minimally invasive and astigmatism-neutral

- **Aberration-free optic⁽¹⁾**

Suitable for all patients regardless of corneal aberration

To preserve depth of field and to enhance contrast sensitivity⁽²⁾

- **360° lens epithelial cell barrier**

To prevent PCO



It is mandatory to follow the instructions for use contained within the packaging

¹⁾ The word 'aberration' as used in this document refers to spherical aberration

²⁾ Compared to aberration correcting and conventional spherical IOLs respectively

DISTRIBUTED BY:

