

Accu-Temp[®] Cautery

Specialty Cauterization

Accu-Temp[®] Cauteries provide a broad selection of tip styles for cauterization in a wide variety of ophthalmic procedures.



For ordering information call:
+44 1865 601 256
or visit www.beaver-visitec.com

 **Beaver Visitec**
Keeping Your Vision in Sight

Accu-Temp[®] Cautery

Specialty Cauterization

High Temperature Cautery




Description	Item #
1205°C (2200°F), 1/2" shaft with fine tip, 10/bx	8442000
1205°C (2200°F), 1/2" shaft with loop tip, 10/bx	8443000
1205°C (2200°F), 2" shaft with loop tip, 10/bx	8444000
1205°C (2200°F), 2" shaft with fine tip, 10/bx	8445000
982°C (1800°F), 1/2" shaft with elongated Vasector fine tip, 10/bx	8441500

Variable Temperature Cautery - Dial Controlled



Description	Item #
454°-649°C (849°-1200°F), 1/2" shaft with fine tip, 10/bx	8440000
813°-1148°C (1495°-2098°F), 1/2" shaft with fine tip, 10/bx	8446000

Low Temperature Cautery



Description	Item #
677°C (1200°F), 1/2" shaft with fine tip, 10/Bx	8441000

Interchangeable Tip Cautery - Reusable Handles, Single-Use Tips



Description	Item #
Variable low temperature handle (requires single-use tips), 10/Bx	8461000
Fixed high temperature handle (requires single-use tips), 10/Bx	8461020
Single-use fine tip, 25/Bx	8461050
Single-use loop tip, 25/Bx	8461060



Beaver-Visitec International Ltd

85c Park Drive, Milton Park,
Abingdon, Oxfordshire OX14 4RY
Tel. : +44 1865 601 256
Fax : +44 1865 595 761
www.beaver-visitec.com