

Visitec® Forceps

Premium Quality Meets Single-Use



- **Precise and consistent tip alignment** for optimal performance with every use
- **Ergonomic handle** providing superior control
- **Safety** — easily identifiable as single-use
- Developed for **MICS**



Capsulorhexis Forceps

Tying Forceps

Bonn Forceps
(Tissue Grasping)

Capsulorhexis Forceps



Finely sharpened teeth enable **precise lens capsule cutting** (581488 only)


Wide platforms allow **excellent lens capsule** grasping and control of the capsulorhexis tear

Perfect alignment

Designed to fit through a **1.8 mm MICS incision**

Sharp and Grasping* tip options

Bonn Forceps (Tissue Grasping)




Excellent tissue grasping and grip strength

Blunt teeth evoke **minimal tissue trauma**

Precise platform alignment to grip sutures

Tying Forceps



Rounded platform edges ensure **sutures are not cut**


Fine tips for **accuracy**

Perfect platform alignment to grip sutures



Our mission is to enable your transition from reusable to single-use instruments, in order to benefit from the many advantages. BVI uses innovative production technology to supply high quality single-use instruments at affordable prices.

- **Improved patient safety**
- **Consistent performance**
- **Time saving and greater convenience**
- **Optimisation of cost and budget control**




Sharp tips




Grasping tips

* Grasping tips should only be used after the capsule has been penetrated by a cystotome or femto laser.



Providing a full range of single-use instrumentation for all your surgical needs

Visitec® single-use instruments are available for insertion into BVI kits.



CustomEyes® Kit



Visitec® Single-use Instrument Tray

For a complete listing of the Visitec® single-use instrument range, please refer to our catalogue at www.beaver-visitec.com

Ordering Information*

Description	Item #	Quantity
Capsulorhexis Forceps with Grasping Tips	581486	10
Capsulorhexis Forceps with Sharp Tips	581488	10
Bonn Forceps (Tissue Grasping)	581489	10
Tying Forceps	581490	10

* Visitec® single-use instruments are available for insertion into BVI kits.



Beaver-Visitec International, Sales Limited

85c Park Drive, Milton Park,
Abingdon, Oxfordshire, OX14 4RY, UK

Tel: +44 (0)1865 601256 Fax: +44 (0)1865 595761

www.beaver-visitec.com