

## Mid-Term Occlusion – Extend™ Absorbable Synthetic Implants


Extend™ Absorbable Synthetic Implants provide occlusion for up to 3 months for patients experiencing short-term dry eye symptoms. These implants are ideal for the transient dry eye associated with LASIK or other ocular surgeries. Extend implants can also be used for allergy sufferers, contact lens discomfort, climate change or to sustain drug contact time and effectiveness.




### A simple solution

- Simple insertion technique
- No insertion tool needed; just forceps
- Can be used to trial permanent occlusion
- No foreign body sensation
- Made of a dissolvable copolymer of glycolic acid and trimethylene carbonate

## Other Dry Eye products

Punctal Gauge/Micro Dilator	Item #	Quantity
	A14-301	1 instrument

Tying Forceps	Item #	Quantity
	581490	10 pack


## Beaver Visitec Dry Eye Products

### Simple insertion and excellent retention

Discover why BVI is the world leader in punctal occlusion – providing solutions for every stage in your patients' continuum of care.



### Extend™ Synthetic Absorbable Implants


	Extend™ Absorbable Synthetic Implants (0.2mm x 2.0mm)	A13-202	20 implants per box, packaged in 10 individual pairs
	Extend™ Absorbable Synthetic Implants (0.3mm x 2.0mm)	A13-203	
	Extend™ Absorbable Synthetic Implants (0.4mm x 2.0mm)	A13-204	
	Extend™ Absorbable Synthetic Implants (0.5mm x 2.0mm)	A13-205	
	SP Extend™ Absorbable Synthetic Implants (0.4 x 2.0mm)	A13-204	1 pair
	SP Extend™ Absorbable Synthetic Implants (0.5 x 2.0mm)	A13-205	
	Extend™ Absorbable Synthetic Implants (Variety Pack)	A13-200	10 pair box of: 3 pair – .3mm 5 pair – .4mm 2 pair – .5mm

**bvi** Beaver Visitec


Beaver-Visitec International, Sales Limited  
85c Park Drive, Milton Park,  
Abingdon, Oxfordshire, OX14 4RY, UK  
Tel: +44 (0)1865 601256 Fax: +44 (0)1865 595761



Permanent Occlusion  
– Parazol® Punctal  
Occluder



Mid-term Occlusion  
– Extend™ Absorbable  
Synthetic Implants





Partial Occlusion  
– Micro Flow™

**bvi** Beaver Visitec  
Keeping Your Vision in Sight

## Permanent Occlusion – Parasol® Punctal Occluder

Permanent occlusion, though easily reversible, Parasol® is indicated for chronic dry eye conditions. Made of medical grade silicone, one glance at the Parasol® Punctal Occluder and you'll recognize a breakthrough product, offering distinct benefits for you and your dry eye patients.

### Easy to insert



### Easy to remove



### Easy to insert

The Parasol's hollow nose collapses on insertion, eliminating the need for punctal dilation in most cases. It's easier on you and more comfortable for your patients.

### Simple sizing

Two sizes, small and medium, fit virtually all patients. X-small and large sizes are available for unusually sized puncta.

### Excellent retention\*

The Parasol is designed for excellent retention, as demonstrated by an independent clinical study.\*

### Size comparison chart

Traditional plugs	0.2 mm	0.3 mm	0.4 mm	0.5 mm	0.6 mm	0.7 mm	0.8 mm	0.9 mm+
Parasol	X-SMALL	SMALL	MEDIUM	LARGE				

Small and medium sizes should accommodate 95% of all patients. BVI recommends undersizing when in doubt.

\*"Plugs Reduce Dry Eye Symptoms, Improve Vision", Craig McCabe, M.D., Review of Ophthalmology, 2009.

## Permanent Occlusion – Packaging Configurations

Parasol® Sterile Pre-Loaded Punctal Occluders	Item #	Quantity
Occluders are packaged in a convenient 1 pair tray. Each occluder is individually sterile and pre-loaded on a high-quality disposable inserter.		
Parasol Sterile Pre-Loaded X-Small (0.25mm – 0.35mm)	A14-200	1 pair tray
Parasol Sterile Pre-Loaded Small (0.35mm – 0.65mm)	A14-201	
Parasol Sterile Pre-Loaded Medium (0.6mm – 0.85mm)	A14-202	2 occluders total
Parasol Sterile Pre-Loaded Large (0.9mm+)	A14-203	

Parasol® Sterile Quickload™ Punctal Occluders	Item #	Quantity
Combines sterile convenience with bulk economy. Occluders are supplied sterile, pre-mounted on a stand. Reusable inserters sold separately.		
Parasol Sterile QuickLoad Small (0.35mm – 0.65mm)	A14-201QL	8 pair per box
Parasol Sterile QuickLoad Medium (0.6mm – 0.85mm)	A14-202QL	
Parasol Sterile QuickLoad Large (0.9mm+)	A14-203QL	16 occluders total
Parasol Sterile QuickLoad Combo Pack (4 pair Small, 4 pair Medium)	A14-205QL	

## Partial Occlusion – Micro Flow™

Designed to provide a simple, reversible means of partial occlusion, Micro Flow™ Punctal Occluders limit the drainage of tears into the lacrimal duct. Ideal when total occlusion is too much.

Micro Flow™	Item #	Quantity
Individually sterile. Pre-loaded onto a high-quality disposable insertion instrument.		
Sterile pre-load Small (0.4mm – 0.5mm)	A14-401	1 pair tray
Sterile pre-load Medium (0.55 – 0.7mm)	A14-402	
Sterile pre-load Large (0.7mm – 0.85mm)	A14-403	2 occluders total

