

LENS DESCRIPTION

Commercial name : Proyal
Lens type: lens for ocular bandage;
Packaging: 5 lenses in blister
Production system : molding
Raw material: biofilcon a group II
Uv absorber: yes
Water content: 58%
Dk 27
Dk/t 33.8
HD aspherical geometry
BC 8,60
Diameter 14,20
Spherical power + 0,50 D
Color Blue tint
CND: Y210309

Tips and tricks

The proyal geometry can be used on almost totality of the corneas, suggesting a range of corneal curvature from 7.6 mm to 8.20 mm.



Proyal lenses are suitable for corneal epithelial lesions and other corneal suffering in addition to eye drops.

Aimino G. - Francone L. - Gaida C. - I.R.I.O. - Istituto di Ricerca in Oftalmologia - Ivrea
studio effettuato in associazione del prodotto "Ekistil gocce oculari".